

HIGH-PRECISION ION MOBILITY CALIBRATION IMPROVES DIA-PASEF ANALYSIS

Tejas Gandhi¹, Stephanie-Kaspar Schoenefeld², Oliver M. Bernhardt¹, David Schlessinger¹, Sven Brehmer², Gary Krappa², Jan Muntel¹, Lynn Verbeke¹, Lukas Reiter¹

1) Biognosys AG, Wagistrasse 21, 8952 Schlieren (Zurich), Switzerland 2) Bruker Daltonik GmbH, Bremen, Germany

Tejas Gandhi, PhD
Bioinformatics Manager

tejas.gandhi@biognosys.com
www.biognosys.com




INTRODUCTION

Ion mobility (IM) based data independent acquisition has recently become available on timsTOF Pro. Due to its parallel accumulation-serial fragmentation (PASEF) strategy it achieves a high duty cycle and has shown promise in biological applications. We have previously shown that high-precision retention time calibration benefits DIA analysis (Bruderer, 2016). Upon applying this concept

to ion mobility for analyzing dia-PASEF data, we could identify a high number of proteins, competitive to state-of-the-art DIA methods.

Figure 1: High-precision IM in Spectronaut
Ion mobility calibration (left panel) ensures that optimal IM extraction width (red line, right panel) is automatically used for data analysis as shown here with a 400 ng HeLa sample.



CONCLUSIONS

- **A novel workflow for high-precision ion mobility calibration in Spectronaut improves identification (Figure 1) and quantification (Figure 3)**
- **Deep protein coverage achieved even with low sample amounts (Figure 2)**
- **SN 14 improves on its dia-PASEF workflow from SN 13 (Figure 3)**
- **Spectronaut processes timsTOF Pro data significantly faster than acquisition time (Figure 2)**

RESULTS




Figure 2: Triplicate HeLa lysate dia-PASEF runs with varying sample amounts (10ng, 50ng, 400ng) on a 100 min gradient

Top Panel: Average number of proteins identified per injected sample amount with 1% peptide and protein FDR. The library consisting of over 10,000 protein groups was created from 15 fractions using Pulsar search engine with 1% FDR (PSM, precursor and protein)

Bottom Panel: The analysis time per run (in minutes) for the library generation (PASEF) and subsequent, targeted analysis of the 400 ng dia-PASEF runs is significantly lower than the acquisition time (100 min gradient + 20 min overhead)




Figure 3: UPS2 proteins spiked in 50 ng mouse cerebellum background matrix in 5 different steps, 2 replicates each

Top Panel: Average number of proteins identified by different Spectronaut versions. The red line indicates the library size in proteins.

Bottom Panel: Regulation analysis was performed using an unpaired t-test of protein quantities calculated from both MS1 and MS2 (Huang, 2019) in SN 14. For SN 13, the t-test was based on precursor quantities calculated from MS2. Upon sorting all the regulated candidate pairs by p-value, we found 34 UPS2 candidate pairs in the top 100. This is comparable to what we have previously observed in this sample (Muntel, 2019).